


LYMPH CARE CONCEPT

CURALYMPH COMFORT

Special textiles and bras to relieve pressure

LIANE

Top



Sizes:

- S (32-38 A; 32-34 B; 32 C)
- M (38-40 A; 36-40 B; 34-38 C; 34 D)
- L (42-44 A; 40-44 B; 38-42 C; 36-40 D; 36-38 DD)
- XL (44-46 B; 42-46 C; 40-44 D; 38-42 DD)
- 2XL (46 C; 44-46 D; 42-46 DD)

Material:

56% Nylon, 36% Polyester, 8% Spandex

Features | Benefits


- **Wide, padded shoulder straps | Wide integrated bra band**
High level of comfort, relief of shoulders, avoids fabric cutting into the tissue so that lymph flow is not hindered
- **Light shaping effect**
High level of comfort for women with mild or light edema
- **Structured/textured material**
Light massage effect
- **Higher back neckline**
For more stability and coverage and high wearing comfort
- **Integrated cup pockets**
To hold breast forms and shapers in place
- **Hip length top**
Ideal for physical activities or sleeping



Style 44812
Sand/Light Sand

LINDA

Wire-free Bra



3D technology channels

Sizes:

- S (32-38 A; 32-34 B; 32 C)
- M (38-40 A; 36-40 B; 34-38 C; 34 D)
- L (42-44 A; 40-44 B; 38-42 C; 36-40 D; 36-38 DD)
- XL (44-46 B; 42-46 C; 40-44 D; 38-42 DD)
- 2XL (46 C; 44-46 D; 42-46 DD)

Material:

67% Nylon, 33% Spandex

Features | Benefits

- **Innovative bonding technology with integrated cup pockets**
Flat seams prevent cutting into the skin and support lymphatic drainage; pockets hold breast forms and shapers in place
- **Wide shoulder straps**
High level of comfort, relief of shoulders, avoids fabric cutting into the tissue so that lymph flow is not hindered
- **Innovative SensElast® 3D technology on the side panels**
Fabric details under the arm that, when the wearer moves, create a light massaging effect, like gentle lymphatic massage
- **Higher cut sides and back**
More support and therefore a higher level of comfort
- **Soft padded back closure**
More coverage and high wearing comfort



Style 44811
Sand



Innovation at work.

These garments have two innovative technologies from SensElast®:


- A **bonding technology** that joins materials together without the need for seams.
 - A **3D technology** using 3D printing to create raised channels on the fabric, like the ones in the underarm area of these two garments.
- Through these two innovations, a wide range of effects can be achieved to bring ultimate comfort to the wearer.

Measurement Guide – CuraLymph Comfort

Finding the perfect fit

Get the right fit in 3 easy steps. Here’s how to measure for your Amoena CuraLymph bra/garment.

1. Underbust measurement: Measure the circumference around the underbust band while still wearing a bra.
2. For a unilateral surgery, use the half-overbust measurement: Measure from the middle of the chest (sternum) over the fullest part of the remaining breast to the center of the spine. Multiply this measurement by 2 to determine the cup size. (Bilateral mastectomy patients can choose the cup size, and it should be proportionate to the body.)
3. Find the bra size required using the sizing chart below.



These instructions are intended as guidelines only. You may need to choose a slightly different size for the perfect fit.

	Cup	Measurements based on a body height of: 5'-5'6"											
Size	-	XS	S	S	M	M	L	L	XL	XL	2XL	2XL	3XL
Underbust measurement in inches	-	25-26	27-28	29-30	31-32	33-35	36-38	39-41	42-44	45-46	47-49	50-52	53-55
Underbust measurement (bra band)	-	30	32	34	36	38	40	42	44	46	48	50	52
Overbust measurement in inches results in the corresponding cup size	A	30-31	32-33	34-35	36-37	38-39	40-41	42-43	44-45	46-47	48-49	50-51	52-53
	B	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54
	C	32-33	34-35	36-37	38-39	40-41	42-43	44-45	46-47	48-49	50-51	52-53	54-55